Tecnovita
b4 **BH**

GAFAS DE MASAJE SONNO

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902 170 258

Thank you for purchasing our Sonno massage unit.

All of our products undergo rigorous testing to ensure maximum quality.

To use Sonno correctly, first read through the whole manual before using it.



INTRODUCTION

Through various forms of vibration, massage and heat, the Sonno glasses can help to increase the oxygen supply, regulate the ciliary muscle function, improve blood circulation to the eyes, relax your eyes and brain, and stimulate your metabolism. The Sonno glasses use microchip and infrared control technology to operate.

Specific uses:

Relaxing eye massage • Combats facial fatigue (dark circles, wrinkles, bags under the eyes) • Anti-stress • Provides a feeling of wellbeing and mental relaxation.

DESCRIPTION

- The product is a multi-functional hi-tech device for massaging and relaxing the eye and temple area.
- Thanks to an advanced microchip, this product combines different variables (vibration, air pressure, infrared heat, magnetic field and relaxing sounds) with full control via a remote control unit.
- Sonno produces an enveloping vibratory motion that stimulates blood flow by means of various multi-frequencies.
- The Sonno glasses use an infrared based heat emitter that activates the radiated cells, causing the molecules to vibrate, strengthening their combination and improving the regenerative capacity of the skin tissue.
- The product successfully applies the concept of air pressure. The pressure provides a general microchip controlled eye massage. This process reduces eye stress.

- The Sonno glasses incorporate a three dimensional magnetic field. This motion provides a constant effective massage to the main parts of the face and eye muscles, generating micro-flows in the blood vessels around the eyes and stimulating blood circulation.
- Sonno is equipped with an MP3 function (512 Mb memory) and a USB connector.
- The product has a built-in rechargeable battery. This can be charged from an electric socket or a computer.

GENERAL SAFETY INSTRUCTIONS


1. Inspect the Sonno glasses before use to ensure that they are in good condition and that the cables are not damaged. If you notice a defect or missing parts, contact the distributor of the product.
2. Do not use this product if you suffer from glaucoma, retinitis or other serious eye diseases or if you have recently undergone surgery.
3. The device must not be used by children or the disabled without suitable adult supervision.
4. Remove contact lenses and any items worn on your hair before using the appliance.
5. Put the glasses on before switching them on.
6. Do not use massage glasses in damp/humid surroundings and prevent them from getting wet. Do not touch the device with wet hands while it is operating. Under no circumstances must

this product be used in the bath or shower.


7. We recommend keeping your eyes closed for a few seconds and then opening them slowly after a massage session.
8. Do not expose the device to direct sunlight or corrosive gases.
9. Switch the device off when you have finished using it.
10. Do not use the glasses when the battery is being charged.
11. Unplug the appliance before cleaning it.
12. Do not wash the Sonno glasses with water. If necessary, use sterile cotton wool to clean them.
13. Do not expose the device to high temperatures. Store it in a dry place.
14. If you notice that the Sonno glasses are not operating correctly, do not attempt to repair, dismantle or modify them. Unauthorised repairs may cause injury, in addition to breaking the terms of the warranty. Contact your technical assistance service.

PREPARATION BEFORE USE

MAIN UNIT

REMOTE CONTROL

- | | |
|-------------------------|---------------------------------------|
| 1. ON/OFF/Timer | 10. Vibration |
| 2. Up selection key. | 11. Air pressure |
| 3. Left selection key. | 12. Heat |
| 4. Function On/Off mode | 13. Music |
| 5. Right selection key. | 14. Volume |
| 6. Down selection key. | 15. USB / Earphone / Charge connector |
| 7. Duration. | 16. LED indicator |
| 8. Battery indicator | 17. Air release valve |
| 9. Selected function | |

Charging the battery via a wall socket:

1. Open the cover for the USB connector on the remote control unit (15).

2. Plug the remote control into the mains supply using the charger.

Charging the battery via a computer:

1. Open the cover for the USB connector on the remote control unit (15).

2. Plug the remote control into the computer using the USB cable.

Caution: Do not use the glasses when the battery is being charged.

Earphone socket.

1. Open the cover for the USB connector on the remote control unit.

2. Plug the earphones into the remote control.

Download the MP3 songs.

To record MP3 songs, first plug the remote control into the computer using the USB cable. Then press the On/Off/Timer button for 3 seconds to

switch the unit on and transfer the songs. (Songs can only be transferred if the unit is switched on).

USE

To switch the device on, press and hold the On/Off/Timer key (1) for 3 seconds.

The device will start operating and the On/Off indicator will light up. The default settings for the Sonno glasses are as follow: 5 minute session, air pressure and heat massage. The symbols for the activated functions (11,12) will appear on the screen.

To choose another function, use the up (2) or down (6) keys to highlight the desired function and a tick will appear to indicate the selected function (9). To activate the selected function, press the Mode key (4). To deactivate the selected function, press and hold the Mode key (4).

You can also change the mode of the selected function by pressing the left (3) or right (5) selection keys.

Air pressure massage

To choose a desired mode, use the up (2) or down (6) selection keys to select the air pressure massage function, activate it by pressing the Mode button (4) (assuming the air pressure function is not already enabled) and select the desired mode by pressing the left (3) or right (5) selection keys.

Mode 1 / Mode 2 / Mode 3 / Mode 4

Vibration massage

To choose a desired mode, use the up (2) or down (6) selection keys to select the vibration massage function, activate it by pressing the Mode button (4) (assuming the vibration

function is not already enabled) and select the desired mode by pressing the left (3) or right (5) selection keys.
Mode 1 / Mode 2 / Mode 3 / Mode 4 / Mode 5 / Mode 6

Heat

To choose a desired mode, use the up (2) or down (6) selection keys to select the heat function, activate it by pressing the Mode button (4) (assuming the heat function is not already enabled) and select the desired mode by pressing the left (3) or right (5) selection keys.
Temperature low 0 / high 1

Music

To listen to the desired song, use the up (2) or down (6) selection keys to select the music function, activate it by pressing the Mode button (4) (assuming the music function is not already enabled) and select the desired song by pressing the left (3) or right (5) selection keys.

Volume

To adjust the volume, first select the volume function by using the up (2) or down (6) selection keys. To increase

the volume press the right (5) selection key. To lower the volume press the left (3) selection key.

Time

To change the time mode press the On/Off/Timer button (1).

5 minutes – 10 minutes – 15 minutes

To switch the unit off, press and hold the On/Off/Timer button (1).

MAINTENANCE & CLEANING

1. Unplug the appliance and allow it to cool down before cleaning it.
2. Clean it with a slightly dampened sponge or soft cloth. Do not allow the appliance to come into contact with water or any other liquid. Do not immerse the appliance in liquid to clean it. Do not use abrasive cleaners, brushes, gasoline, glass or furniture polishes or solvents to clean it.

TECHNICAL SPECIFICATIONS

MODEL	SONNO YM68
AMBIENT OPERATING TEMPERATURE	5°C - 40°C
AIR PRESSURE	860hPa – 1,060hPa
MAGNETIC FIELD INTENSITY	50mT – 250mT
HEATING TEMPERATURE	30°C – 45°C
POWER SUPPLY	DC: 5V 500mA 6VA
DIMENSIONS: REMOTE CONTROL	14.1 x 5.55 x 3.3cm.
DIMENSIONS: MAIN MASK	20.6 x 10.8 x 8.6cm.
NET WEIGHT	458g

Warranty

This product is guaranteed for 2 years in accordance with that stated in Directive 1999/44/EC of the European Parliament and of the Council 25th May, 1999, regarding certain aspects of the sale of consumer goods and associated guarantees, incorporated into Spanish Law by virtue of Law 23/2003, 10th July; Guarantees on the sale of consumer goods.

To this effect, EXERCYCLE, S.L. VAT nº B-01029032, with business address c/ Zurrupitieta 22, Vitoria, Spain, guarantees this product.

Date of purchase.-	Purchased in.-
Purchaser.-	Dealer stamp.-



Please keep the receipt as proof of purchase.

Should you have any queries please do not hesitate to phone our customer service line:

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BH China Co., Ltd.